

Lungs



IMPORTANT

- 1** In treating the heart area (4th intercostal, left of sternum) NEVER use a frequency other than 5 Hz.
- 2** A number of treatments in this user guide are complex treatments in which several stand-alone treatments have been combined to form one course of treatment. If in such cases more than one of the individual treatments calls for the heart area (4th intercostal, left of sternum) to be treated, ensure that this is done only once.
- 3** Do NOT use coMra therapy on patients who have had heart transplants, who have hearts with artificial valves or who have cardiostimulators.



Treatment times given in this user guide were designed for the home user with the Delta or Palm 980 nm laser to give a conservative starting point. Treatments with the professional 905 nm laser are more intense, and therefore these treatment times may be shortened, and also the space between treatments extended further, as determined by feeling and experience.



In treating the heart area NEVER use a frequency other than 5 Hz.


Do NOT use on patients who have had heart transplants, who have hearts with artificial valves or who have cardiostimulators.

SECTION 7

PULMONOLOGY 1

24 mins

Treatment time.



Using:
Delta Medical Terminal or
coMra Palm.


1 x daily
for 15 x days

- PNEUMONIA
- PLEURISY
- BRONCHITIS

TREATMENT POINT	⌚ MINS	⏏ HZ
1 4 th intercostal, left of sternum (only once per day)	2	5
2 Subclavian fossae	2 per point	
3 Supraclavicular fossae	1 per point	
4 3 points along both sides of vertebral column, opposite scapulae	1 per point, per frequency	5 then 50
5 2 points below scapulae	2 per point	5



Treatments with the professional 905 nm laser are more intense, and therefore these treatment times may be shortened, and also the space between treatments extended further, as determined by feeling and experience.




SECTION 7

PULMONOLOGY 2

23 mins

Treatment time.



Using:
Delta Medical Terminal or
coMra Palm.

1 x daily

for 15 x days


Allow 3 weeks rest
and repeat cycle until
improvement is noticed.

• ASTHMA

TREATMENT POINT	⌚ MINS	⚡ HZ
1 Manubrium sterni	2	5
2 Sternum	1	
3 Trachea in thyroid region	1 per frequency	5 then 50
4 Subclavian fossae	1 per point	5
5 Supraclavicular fossae		
6 Right subcostal	1	50
7 Left subcostal		
8 Epigastrium		
9 3 points across base of neck	1 per point	1000
10 3 points along both sides of vertebral column, opposite scapulae		50
11 Point between thumb and index finger on back of both hands	1 per hand	1000



Treatments with the professional 905 nm laser are more intense, and therefore these treatment times may be shortened, and also the space between treatments extended further, as determined by feeling and experience.





In treating the heart area NEVER use a frequency other than 5 Hz.


Do NOT use on patients who have had heart transplants, who have hearts with artificial valves or who have cardiostimulators.

SECTION 7

PULMONOLOGY 3

48 mins

Treatment time.



Using:
Delta Medical Terminal or
coMra Palm.

1 x daily
for 21 x days

Allow 2 weeks rest and
repeat cycle until disease is
completely cleared.

• TUBERCULOSIS

MORNING Apply UNIVERSAL 3 (BLOOD) every day

EVENING	TREATMENT POINT	⌚ MINS	⚡ HZ
1	4 th intercostal, left of sternum (only once per day)	2	5
2	Subclavian fossae	2 per point	
3	Supraclavicular fossae		
4	4 points along both sides of vertebral column	1 per point, per frequency	5 then 50
5	3 points over each lung	2 per point	5



Treatments with the professional 905 nm laser are more intense, and therefore these treatment times may be shortened, and also the space between treatments extended further, as determined by feeling and experience.

