



# DeltaSeries Battery Operating Instructions

*Please Note: Due to the Samsung Note 7 devices causing fires manufacturers have had to change the way batteries operate. These detailed instructions reflect these changes.*

1. The battery is shipped with a storage charge only. **You must charge the battery before use.**
2. To charge the battery; plug the small end of the battery charging cable into the battery charging port on the battery. Plug the large end of the battery charging cable into the wall plug adapter. Plug the wall plug adapter into a power outlet.
3. Once plugged into the power outlet the indicating light will turn red within 1-5 seconds. The indicating light will remain red as long as the battery is plugged in to charge. (On some models the red light will blink, on others it will be steady.)
4. Allow the battery to charge for no less than 12 hours continuously before use. After you have completed the initial 12 hour charge the battery is ready for operation. **Failure to do this initial charge may reduce the battery's life.**
5. To operate your Delta device using the battery; plug the Delta device USB power cord into the battery power port for the Delta device. Once plugged in the indicating light will turn blue.
6. You will have approximately 30 seconds to operate your Delta device once plugged into the battery. If you do not start the treatment program on your Delta device within 30 seconds the battery will automatically shut off and the (blue) indicating light will go out.
7. If this happens you will need to unplug the Delta device USB power cord from the battery, wait 5-7 seconds and then plug it back in to operate the Delta. If you try to plug the Delta device USB power cord into the battery before the 5-7 second pause, the indicating light will not light up blue and the battery will not power up the Delta device.

